

Validations and Certifications

nPLEX Food Total kit





Catalogue number: TAGN01001

VALIDATIONS

nPLEX Food Total has been validated for the detection of *Escherichia coli*, *Listeria monocytogenes*, *Salmonella spp.* and *Staphylococcus aureus* in several artificially contaminated food samples, such as:

CATEGORY	CONDITION	FOOD SAMPLE
Poultry	Raw	Whole chicken, Turkey breast, Chicken fillet, Turkey ribs.
	Heat-processed	Roasted turkey breast, Roasted turkey legs, Whole roasted chicken, Roasted chicken breast
	Frozen	Frozen chicken fillets, Whole frozen turkey
Meat	Raw	Ground beef, Minced beef, Beef tenderloin, Pork tenderloin, Ground pork
	Heat-processed	Roasted meat and rice, Roasted chicken and mashed potatoes
	Frozen	Frozen pork meat
	Raw	Fruit Mix; Cabbage, broccoli and cauliflower, Mixed salad
Fruits and vegetables	Heat-processed	Pasteurized apple juice
	Frozen	Frozen broad beans, Frozen corn, Frozen broccoli
	Fermented	Fermented cabbage
	Cured / salted	Pickles
	Other	Mushrooms
Fish and Sea-Products	Raw	Fresh salmon
	Heat-processed	Seafood assortment
	Frozen	Frozen tilapia
	Smoked	Smoked salmon
Dairy	Raw	Raw milk
	Heat-processed	Pasteurized milk
	Frozen	Chocolate ice-cream
	Fermented	Strawberry yoghurt
	Dried	Powdered milk 26% fat
Chocolate and Bakery	Low humidity	Chocolate 85% Cocoa
	Dried	Cake mix
Other Products	Spices	Oregano
	Sauces	Mole sauce
	Eggs	Fresh eggs

CERTIFICATIONS

SYMBOL	MEANING
	In Vitro Diagnostics medical device
	This product fulfills the requirements of the European Directive 98/79 EC for <i>in vitro</i> diagnostic medical devices
	Authorized Representative in the European Community
	The technique upon which this product was designed has been accredited by the American Association for Laboratory Accreditation